


Cosmetics

Comfort

Control



Radix®



Genux® OA



Wilmer®



Ambroise Products are designed according to the 3C philosophy. This design philosophy states that a good prosthesis or orthosis should comply to at least three major characteristics: good cosmetics, maximal comfort and optimal control. The UTX® orthoses fulfil these requirements, as do the Radix® and Genux® OA. The Wilmer® line of orthoses is designed according to this philosophy as well. With this the user of Ambroise products is placed at the center of focus.

Ambroise

that fits your body dat past bij je lichaam


The dynamic UTX® orthoses



Ambroise

that fits your body dat past bij je lichaam


UTX[®]-Swing



The UTX[®] is the line of dynamic KAFOs that enables the user to walk naturally. The UTX[®] supports when needed during gait in cases of muscular insufficiency.

Indication:
Independent of pathology. The amount of rest function in hip and knee extensors determines the most appropriate type of UTX[®]. See also the UTX[®] selection protocol.


UTX[®]-FS



The UTX[®]-FS support the knee in the coronal plane in cases of valgus or varus instability. The UTX[®] FS offers optimal stability during gait and reduces pain in the knee.

Indication:
Equal to the UTX[®] SWING, but added with required knee stability in the coronal plane.


UTX[®]-Post



The UTX[®]- Post supports the knee in the sagittale plane in cases of hyperextension. The UTX[®] POST prevents increase of the hyperextension and corrects the leg to neutral alignment.

Indication:
Equal to the UTX[®] SWING, but added with required sagittal control in cases of Genu Recurvatum.

UTX[®]-DZ



The UTX[®]-DZ supports the ankle in the coronal plane and offers support in dorsal flexion. This makes the UTX[®] a line of KAFOs with a suitable solution that fits the body for a wide variety of pathologies.

Indication:
Equal to UTX[®] SWING, but added with required stability in the sagittal and frontal plane around the ankle.

Standard colour choices: Neutral coloured frame with either black or beige straps and pelottes.